

SM-ZQ12

Classic 12-Ch ECG Machine




**S**  
**Summa**  
**Medical**



Features

- 12-lead synchronous acquisition and display
- 9.7-inch tilted touch screen with grid background
- High-performance chip module makes ECG signal collection more accurate
- Sampling rate: 32000Hz
- A/D shift: 24 bits
- 0.01-350Hz ultra-wide frequency response could acquire signals of patients from pediatric to elder
- With waveform freezing and waveform review functions, clinicians can choose any waveform for analysis and printing
- With report preview function, clinicians can ensure signal quality
- Support reanalysis based on modified patient information
- Support USB/ RJ45/ SD card ports for peripherals and export of data
- Support output of multiple formats

Network and Connectivity



Innovative Ergonomics



Report Function

