

Features

- 5-inch tilted screen with grid background
- 12-lead synchronous acquisition and display
- Sampling rate: 32000Hz
- 0.01-350Hz ultra-wide frequency response could acquire signals of patients from pediatric to elder
- With report management function, it supports report preview, query, editing, printing and transmitting
- Support paperless records and storage of electronic cases, the machine stores more than 1000 cases
- Support output of multiple formats
- Support wired and wireless way to connect Internet
- Optional lifting trolley



Waterproof and pressure-proof suitcase
(optional)

SM-ZQ3

3-Ch ECG Machine

