

# SM-OT001 Economic Manual Operating Table

This product is intended for use in hospital operating rooms for surgical procedures involving the head, neck, thoracic and abdominal cavities, perineum, and limbs, as well as for gynecological, ophthalmic, ENT, and orthopedic surgeries.



Item	Parameter
Tabletop width	(2100×500)±20mm
Tabletop height	Min750 (±20)mm
Height adjustment	250mm
Tabletop forward tilt	Forward tilt≥30° Backward tilt≥15°
Lateral tilt	≥18°
Head plate adjustment	Up≥45° Down≥90°
Back plate adjustment	Up≥70° Down≥18°
Waist plate adjustment	≥120 mm
Leg plate adjustment	≥90°

- The angle adjustment of the backrest and leg plates of this general surgical table is achieved through a hand-cranked mechanism using worm gears and rack-and-pinion transmission.
- The tabletop and covers can be made of high-quality cold-rolled carbon steel with powder coating or premium stainless steel, depending on user requirements.
- Hydraulic lifting is driven by an oil pump, and all required patient position adjustments during surgery are operated from both sides of the tabletop.

