

## Model :SM-W809 Steel Wheelchair

- ★ Chromed Steel wheelchair
- ★ Seat Cushion :Oxford/PVC
- ★ Fixed PVC Armrest and Fixed Footrest
- ★ 8 Inch PVC Wheels
- ★ 24" Spoke Rear Wheel with Solid Tire.
- ★ Aluminum Footplate
- ★ Weight Capacity :100kg
- ★ N.W:15KG
- ★ G.W:17KG
- ★ Packing Size:  
91.5\*21\*89cm.1pcs/ctn

PS:  
Manual measurement  
of dimensions.tolerance  
range  $\pm 1$ cm

