

MODEL SM-W867

ALUMINIUM FRAME LIGHTWEIGHT AND FOLDABLE

Self-locking brakes/effortless pushing/multi-colour options



Aluminium frame, double cross bar, Drop back handle brake (with self-locking function, Fixed Inclined, Flip up Footrest seat back pad (optional breathable mesh), 6" solid front wheel, 20" PU steel wire rear wheel, rear detachable push handle, plastic footrest



MAIN PARAMETERS

Total length: 96cm
Total width :62cm
Total height: 88cm
Package size: 79*27*72cm
Folded width: 27cm

Seat width: 46cm
Seat depth: 40cm
Seat height: 49cm
Height of backrest: 42cm
Handrail height: 19.5cm

Front wheel diameter: 15cm
Rear wheel diameter: 51cm
Maximum load capacity: 100kg
Net weight: 11kg