

## BLACK STEEL HEAVY DUTY WHEELCHAIR

Designed for obese patients with larger body weights

SM-W909



### Product Details



Flip up armrest



Detachable footrest  
Black plastic footplate



8 inch solid front wheels



24" MAG Rear Wheels with Solid Tyre.



Weight capacity 150kg

### Product Dimensions

Total length	110cm
Overall width	70.5cm
Overall height	93cm
Seat width	55cm
Seat depth	40cm
Seat height	52cm
Front wheel	8 inch solid wheels
Rear wheel	24" MAG Rear Wheels with Solid Tyre.
Folding width	28cm
Max capacity	150kg

