



- ① The operation interface adopts heat-sensitive touch buttons, which is sensitive to operation and has a long service life.
- ② Height is adjusted by hydraulic system, ergonomic height and convenient for nursing. HI-LOW: 960-840mm; Trend./Reverse Trend.: 0-10°, effectively prevent the baby from milk regurgitation problem etc.
- ③ The handrail is easy to operate and adjustable. When the handrail is vertical to 90°, it can be used as the support of mosquito net
- ④ Humanized design:
Music playback function, using pluggable SD card, user can change music at will Weighing function: measure the baby's weight at any time Indoor temperature measuring
- ⑤ Net weight: 27KG



- ① 操作界面采用热感应触摸按键，操作灵敏使用寿命长
- ② 通过液压系统调节车体高低，符合人体工程学高度方便护理。高低调节范围：960-840mm；具有倾斜功能，角度调节范围 0-10°，有效的防止宝宝吐奶问题
- ③ 扶手操作方便，可调节。当扶手垂直到 90° 时，可作为蚊帐的支架使用
- ④ 人性化设计
音乐播放功能：采用可插拔 SD 卡，用户可以随心更换音乐
称重功能：可随时测量婴儿的体重
可实时测量室内温度
- ⑤ 净重：27KG