

# Constant-Temperature Incubator



## Features:

- \* The inner chamber is made of high-quality brushed stainless steel with a semi-circular four-corner design for easy cleaning.
- \* All inner surfaces are wrapped with aluminum foil, and the gaps between the inner liner and partitions are filled with high-quality insulation cotton, providing excellent heating and heat retention for more uniform temperature distribution.
- \* Equipped with leakage current and over-voltage protection, over-temperature alarm, and a temperature control switch that cuts off power when overheating occurs to prevent damage to the instrument.

## Technical Parameters:

| Model                   | BJPX-H30L                                           |
|-------------------------|-----------------------------------------------------|
| Capacity                | 30L                                                 |
| Temperature Sensor      | PT100                                               |
| Temperature Range       | RT+5~85°C                                           |
| Temperature Fluctuation | ±0.2°C (@37°C)                                      |
| Temperature Uniformity  | ±0.5°C (@37°C)                                      |
| Temperature Precision   | 0.1°C                                               |
| Timing Range            | 0~9999min/h                                         |
| Display                 | LCD                                                 |
| Power Supply            | AC 220V, 50/60Hz(Standard); 110V, 50/60Hz(Optional) |
| Consumption             | 150W                                                |
| Shelves                 | 2pcs                                                |
| Standard Accessories    | Shelf                                               |
| Optional Accessories    | Shelf, footmaster caster                            |
| Internal Size(W*D*H)    | 310*310*310mm                                       |
| External Size(W*D*H)    | 432*538*570mm                                       |
| Package Size(W*D*H)     | 550*670*700mm                                       |
| Net Weight              | 28kg                                                |
| Gross Weight            | 40kg                                                |