

Children Cartoon Stethoscope

儿童型卡通听诊器



Champaign gold



Dumb silver



Dumb gun color

Patent protection: ZL202222660847. 5

Patent protection: ZL202121114900. 0

Patent protection: ZL202130418193. 3



Royalblue



Purple



Yellow



Red



Dark Blue

SM-SC001

Children Cartoon Stethoscope

Material(材料)

- Chestpiece: ABS / ABS(Vacuum Plating)+ Zinc Alloy
- Binaural: Stainless steel
- Y tube: PVC

- 听头: ABS / ABS(真空电镀)+锌合金
- 耳挂: 不锈钢
- Y管: PVC

Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
260x130x38	675x275x415	50	14/13